READ CAREFULLY BEFORE USE

Light-cured Dental restorative material

BEAUTIFIL II LS is a light-cured (400-500 nm), radiopaque restorative material indicated for anterior and posterior restorations including occlusal surfaces.

COMPOSITION

Glass powder, Urethane diacrylate, Bis-MPEPP, Bis-GMA, TEGDMA, Polymerization initiator, Pigments and others

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➢ Repair of restorations and prostheses

PRECAUTIONS

- Do not use this product on patients with known allergies to this material and/or methacrylate monomer.
- Operators with known allergies to this material and/or methacrylate monomer should exercise caution while using this product.
- If any inflammation or other allergic reactions occur on either patient or operator, immediately discontinue use and seek medical advice.
- Use medical (dental) gloves and protective glasses to avoid sensitization to this product. In case of accidental contact with oral soft tissue or skin, immediately blot with alcohol moistened cotton ball, and rinse with plenty of water. Medical (dental) gloves provide some protection against direct contact with this material; however, some monomers can penetrate through the gloves very quickly to reach the skin. Therefore, should this material contact even the gloves, immediately remove them and wash hands with plenty of water.
- In case of contact with eyes, immediately flush the eyes with plenty of water and seek medical advice.
- Do not re-use tips.
- This product is intended for use by dental professionals only.

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A-1. Cleaning of tooth surface

Thoroughly clean the tooth surface to remove plaque.

A-2. Shade selection

Select appropriate shade using a shade guide while the tooth is still moist.

A-3. Cavity preparation

Prepare the cavity following the standard method.

A-4. Isolation

Use a rubber dam etc to isolate the treatment area from contamination.

A-5. Placement of matrix

Place a matrix band in case the cavity involves proximal area. Follow the matrix manufacturer's instructions for use.

A-6. Pulp protection

In case of pulp exposure or deep cavity with close proximity to the pulp, cap the pulp with calcium hydroxide etc.

A-7. Adhesive application

Apply dentin adhesive following the manufacturer's instructions for use. Either two step adhesive such as FL-BOND II or one step adhesive such as BeautiBond can be used.

A-8. Flowable application

To achieve optimal adaptation of the prepared tooth, apply a flowable restorative material such as BEAUTIFIL Flow or BEAUTIFIL Flow Plus as a liner.

A-9. Filling

Dispense adequate amount of BEAUTIFIL II LS onto a paper pad from the syringe. Apply the dispensed material into the cavity using a suitable instrument and create the shape desired. When using a tip, insert the tip into the dispenser to apply the material directly into the cavity. The tip material may also be dispensed onto a paper pad to be applied into the cavity with a suitable instrument.

➢ When storing these materials in a refrigerator, allow them to reach room temperature before use.

➢ The paste dispensed onto the paper pad should be shielded from ambient light with light-shielding cover to avoid premature polymerization of the material.

➢ Tightly cap the syringe after each use.

➢ When applying the material directly into the cavity, sterilize the dispenser before and after use.

➢ When the cavity is deep, incremental filling and light-curing of the material in layers, not exceeding 2 mm each, is recommended.

A-10. Light-cure

Light-cure with a dental light-curing unit. In cases where the restored area is large, light-cure in segments.

Light-curing unit	Halogen	LED
Light-curing time	20 sec.	10 sec.
Wavelength	400-500 nm	440-490 nm
Light intensity	≥ 500 mW/cm²	≥ 1000 mW/cm²

➢ Use a dental light-curing unit having specifications equivalent to those described above. When using the one having performance scale other than these recommendations, follow the instructions of each unit light-cure BEAUTIFIL II LS properly.

➢ Operators and patients must use protective glasses or light shielding plate to avoid direct exposure to the curing light.

➢ Disinfect the light-curing unit with alcohol after each use to avoid cross-contamination.

A-11. Finishing & polishing

After light-curing, contour, finish and polish the restoration following the standard method.

B. Repair of restorations and prostheses

B-1. Cleaning of the restoration to be repaired

Thoroughly clean the restoration to remove plaque.

B-2. Shade selection

Select appropriate shade using a shade guide.

B-3. Roughing of the restoration

Roughen the surface to be repaired with a dental abrasive (diamond burs etc). Then, rinse with water and dry.

B-4. Isolation

Use a rubber dam etc to isolate the treatment area.

B-5. Pretreatment of the restoration

Pretreat the material to be repaired (composite resin, porcelain, metal or natural tooth) to enhance mechanical and chemical retention.

B-6. Filling, light-cure, finishing and polishing

Follow the procedures from A-9 to A-11 above.

NOTE ON USE

- Do not use eugenol containing materials since they might inhibit polymerization or adhesion of the material.
- Follow the instructions for use of each dental material, instrument or equipment to be used in conjunction with this product.
- Do not use this product for any purposes other than specifically outlined under INDICATIONS in these instructions for use.
- Use this product within the expiration date indicated on the package. (Example □ YYYY-MM-DD → Year-Month-Date of the expiration date)

STORAGE

Store in a dark place at room temperature (1-30°C/34-86°F). Avoid high temperatures and high humidity. Keep away from direct sunlight and any source of ignition.

CAUTION: US Federal law restricts this device to sale by or on the order of a dental professional.

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