

## BEAUTICEM SA



PRG  
Technology  
Giomer

- Self-etching, self-adhesive resin cement
- S-PRG (Surface Pre-reacted Glass Ionomer) – ability to recharge fluoride
- No primer required (except porcelain)
- High bond values across all substrates
- Auto-mixing & dual cure
- High levels of radiopacity
- Low film thickness (12µm)

## BEAUTICEM SA

PN

Product

**3213 . . . . . BeautiCem SA Clear**

5ml paste, 10 mixer tips & step card

**3214 . . . . . BeautiCem SA Ivory**

5ml paste, 10 mixer tips & step card