

# CX-Plus Capsule

Glaslonomer Luting Cement



# CX-Plus Capsule

Glaslonomer Luting Cement

CX-Plus Capsules are designed to save time and provide superior handling. Indicated for cementing crowns, inlays, long span bridges, and orthodontic appliances, CX-Plus has a high compressive strength, low film thickness, and is radiopaque. It is resistant to water contamination and can exist in a moist environment.

- High compressive strength at 149 MPa compared to other leading competitors.
- Low film thickness at 15 microns optimizes marginal integrity
- Radiopacity at 1.55mm
- Low levels of acid erosion at 0.02mm
- High fluoride release
- 0.53g per capsule
- Work/Set time 4 minutes 10 seconds from start of mixing
- 3 year shelf life

"I have used CX-Plus Luting Cement for stainless steel crowns and zirconium crowns for children. This cement is one of my favorites, it saves time and is reliable. Shofu's CX-Plus Luting Cement is great for pediatric crowns and works well, even for difficult moisture-control situations."

- Fred Margolis, DDS, FICD, FACD, MALD  
Buffalo Grove, IL

Visit [www.shofu.com](http://www.shofu.com) or call **800.827.4638**



One: Before



Two: After



Photos courtesy of Fred Margolis, DDS  
NuSmile Pediatric Crowns

CASE STUDY



#### Order Information

Part Number 3277  
50 Pre-mixed Capsules  
Cx-Plus still available in  
powder/liquid form  
Part Number 1166



Shofu Dental Corporation  
1225 Stone Drive  
San Marcos, CA  
92078-4059, USA

SN196-113

Glaslonomer Luting Cement