

## CX-PLUS GLASIONOMER (TYPE I LUTING CEMENT)



- NOW IN CAPSULES!
- Indicated for cementing crowns, inlays, long-span bridges & orthodontic appliances
- Low film thickness (15µm) with high compressive strength (149 MPa)
- Bonds to both dentin & enamel to provide an excellent marginal seal to overcome post-operative sensitivity
- High fluoride release & high level of radiopacity
- Easy mixing, handling & resistant to water contamination at the initial stage
- 3 year shelf life

### CX-PLUS GLASIONOMER

#### PN

#### Product

**1166** ..... **CX-Plus Cement Kit**

35g powder, 17ml liquid, mixing pad, spatula, dispenser & instructions

**1167** ..... **CX-Plus Powder 35g**

**1168** ..... **CX-Plus Liquid 17ml**

**1169** ..... **CX-Plus Triple Kit**

3 (35g) powder, 3 (17ml) liquid, mixing pad, spatula, dispenser & instructions

**3277** ..... **CX-Plus Capsules**

50 (0.53g) capsules & instructions