

Super-Snap **X-TREME**TM

FINISHING & POLISHING
ULTRA-GLOSS PERFORMANCE KIT

BONUS INSIDE!

Beautiful Flow Plus
& STORAGE CASE

SHOFU



SN242-0215



Super-Snap X-Treme
Technique Video

Visit www.shofu.com or call **800.827.4638**

Shofu Dental Corporation • San Marcos, CA


Super-Snap X-TREME

FINISHING & POLISHING
ULTRA-GLOSS PERFORMANCE KIT


For X-Treme Ultra-Gloss Results

The Ultra-Gloss Performance Kit incorporates the essential tools to complete quick, easy and highly aesthetic restorations. A system guaranteed to provide patients with the wet, glossy look they desire. Kit contents include 225 disks and more!

- **4-step, color coded disk system for contouring, finishing and polishing**
 - Integrating Super-Snap coarse/medium disks with Super-Snap X-Treme fine/superfine disks
 - X-Treme disks offer greater thickness for increased durability, tactile feel and flexibility without warping
 - New 3D semispherical grit reduces clogging, secondary scratches and achieves X-Treme gloss results
 - Gouging eliminated – no metal center
 - Two sizes - Four grits (aluminum oxide and silicon carbide)
- **DirectDia Paste - a 20% diamond polishing paste for high luster polishing**
 - Easily distinguishable green color and a pleasant lime flavor
- **Super-Snap Super Buff Disks – synthetic velvet clothed disk**
 - Easily retains paste while polishing
 - Very light touch needed to produce a lustrous shine



400x 200µm
3D semispherical coating
Super-Snap X-Treme Red Disk



*Test results on file. All trademarks belong to their manufacturers.

Super-Snap X-Treme produced 41% greater gloss than Sof-Lex Disks and 29% greater gloss than Sof-Lex Spirals



NEW! Super-Snap X-Treme Ultra-Gloss Performance Kits

Standard Ultra-Gloss
Performance Kit 12mm
(PN 0500XS)

Contents (225 pieces):

Super-Snap Disks:

L506 50 pcs

L528 50 pcs

Super-Snap X-Treme Disks:

L511 50 pcs

L512 50 pcs

Plus:

L523 Super-Snap Buff Disks, 25 pcs

0439 Super-Snap CA Mandrel 4 pcs

0558 1 DirectDia Polishing Paste, 3g

BONUS:

2002 1 Beautifil Flow Plus

F00, A2, 2.2g syringe

6 compartment storage case

Mini Ultra-Gloss
Performance Kit 8mm
(PN 0500XM)

Contents (225 pieces):

Super-Snap Disks:

L507 50 pcs

L522 50 pcs

Super-Snap X-Treme Disks:

L531 50 pcs

L532 50 pcs

Plus:

L524 Super-Snap Buff Disks, 25 pcs

0439 Super-Snap CA Mandrel 4 pcs

0558 1 DirectDia Polishing Paste, 3g

BONUS:

2015 1 Beautifil Flow Plus

F03, A2, 2.2g syringe

6 compartment case storage

SHAPE CATEGORY	Polishing		Super-Polishing	
Super-Snap X-Treme Shape				
Color/Grit	Green (fine)		Red (super fine)	
Refills (50 pk)	PN	L511 L531	L512	L532

SHAPE CATEGORY	Contouring		Finishing			
Super-Snap Shape						
Color/Grit	Black (coarse)		Violet (medium)			
Refills (50 pk)	PN	L506 L507	L508*	L528	L509*	L522

*Not included in the Super-Snap X-Treme Ultra-Gloss Performance Kits