











NEW! ZiLMaster POINTS



- ZiLMaster is a three-step system for contouring, finishing, pre-polishing & polishing
- Diamond impregnated silicone polishers that come in three different grits
- Coarse is suitable for abrasive prepolishing & finishing of areas or zones that still need minor material reduction
- Medium is recommended for contouring & smoothing anatomical details & preparing surfaces for high-gloss polishing
- Fine is designed for efficient high-gloss polishing without the additional use of polishing pastes
- **Use ZiLMaster for:** zirconia & lithium disilicate

| SHAPE CATEGORY | Bullet | Knife | Cup | Mini-Pt. | KN7 | WH6 | Unm. Wheel | Kits |
|----------------------------|---|---|---|---|---|---|---|-------------------|
| SHAPE |  |  |  |  |  |  |  | See Assort. Below |
| ISO# | 040 | 090 | 065 | 030 | 125 | 120 | 170 | |
| COARSE (green band) | | | | | | | | |
| HP (3 PK) | PN 0640 | - | - | 0641 | 0642 | 0643 | 0645 (1 per pack) | 0644 |
| CA (3 PK) | PN 0658 | 0662 | 0663 | 0661 | - | - | | 0664 |
| MEDIUM (blue band) | | | | | | | | |
| HP (3 PK) | PN 0646 | - | - | 0647 | 0648 | 0649 | 0651 (1 per pack) | 0650 |
| CA (3 PK) | PN 0665 | 0659 | 0668 | 0666 | - | - | | 0669 |
| FINE (yellow band) | | | | | | | | |
| HP (3 PK) | PN 0652 | - | - | 0653 | 0654 | 0655 | 0657 (1 per pack) | 0656 |
| CA (3 PK) | PN 0670 | 0672 | 0673 | 0671 | - | - | | 0674 |

ZiLMaster POINTS

PN

Product

0644 **Coarse Assortment Package - HP Shank**
 Bullet, Mini-Point, KN7, WH6

0650 **Medium Assortment Package - HP Shank**
 Bullet, Mini-Point, KN7, WH6

0656 **Fine Assortment Package - HP Shank**
 Bullet, Mini-Point, KN7, WH6

ABRASIVES

NEW! ZILMASTER POINTS

ZILMASTER POINTS

| PN | Product |
|-----------|--|
| 0664 | Coarse Assortment Package - CA Shank Bullet, Knife, Cup, Mini-Point |
| 0669 | Medium Assortment Package - CA Shank Bullet, Knife, Cup, Mini-Point |
| 0674 | Fine Assortment Package - CA Shank Bullet, Knife, Cup, Mini-Point |
| 0680 | Adjustment Kit - HP Shank |
| 0681 | Adjustment Kit - CA Shank |
| See chart | Coarse CA & HP Shank (3 per pack) |
| See chart | Coarse Unmounted Wheel (1 per pack) |
| See chart | Medium CA & HP Shank (3 per pack) |
| See chart | Medium Unmounted Wheel (1 per pack) |
| See chart | Fine CA & HP Shank (3 per pack) |
| See chart | Fine Unmounted Wheel (1 per pack) |